

HemoCleanse Technologies, LLC

3601 Sagamore Parkway North, Suite B

Lafayette, Indiana 47904-5033 USA

Phone: (765) 742-4813

FAX: (765) 742-4823

Website: www.hemocleanse.com

Ordering: dcarr@hemocleanse.com

Price List

IJ Catheter Ultrasound Model, Version 2400

September, 2017

Model Purchase Cost: **\$10,000**

Standard set of replacement parts: 40 silicone rubber subcutaneous layers, 80 veins, 10 arteries and 4 O-rings for each vein and artery.

Standard set of spare parts:

O-ring for vascular insert

Heart Bulb with Tube

Neck – Reservoir Tube

Tunnel Site cover, plain (2)

Replacement Parts and Optional Supplies:

Part Name	Quantity	Price
Subcutaneous Layer, Red Silicone Rubber	1	\$6.00
Subcutaneous Layer, Clear Silicone	1	\$16.00
Skin Assembly	1	\$60.00
Vein	1	\$1.40
Artery	1	\$2.80
O-ring for Vein, # 18	100	\$13.00
O-ring for Artery, # 12	100	\$9.50
O-ring for Vascular Insert, # 241	1	\$4.20
Tunnel Site Cover, Velcro, plain	1	\$2.60
Tunnel Site Cover, Velcro with slot	1	\$5.20
Heart Bulb with tube	1	\$42.00
Neck - Reservoir Tube	1	\$21.00
Water Reservoir	1	\$125.00
Water Reservoir Holder	1	\$180.00
Site-Rite Needle Guide Kit, 18G	1	\$30.00
Site-Rite Needle Guide Kit, 21G	1	\$30.00
Ultrasound Gel, 250 mL bottle	1	\$6.00
Syringe, 10 mL	1	\$0.20

Notes: Skin and veins are consumed during use. Each skin and set of veins should be good for at least ten catheter insertions on each side.

O-rings are reusable; Arteries are usually not consumed during use.

Terms:

- ❖ Shipping and handling cost is billed in addition to sale or rental.
- ❖ Domestic Shipments: 50% with purchase order and 50% net 30 days.
- ❖ International Shipments:
 - Items will be shipped after payment for the goods and shipping & handling is received in full.
 - Customer pays for all customs fees, duties and taxes. Preferred shipping terms are DAT: Delivery At Terminal. If import fees or taxes are well defined and low, DAP or DDP are acceptable.
 - Preferred shipping company is Federal Express.

Rental:

Rental Price: \$1,500 per model, plus shipping and handling, for one week of use.

Payment Terms: Domestic: Net 30 days. International: Full payment in advance.